


## ModDetect<sup>™</sup> Panels

# Detect nucleic acid modifications independent of sequence or location

The first-of-its-kind specialty reagent panels designed to streamline development of oligonucleotide therapeutics, mRNA vaccines, and more.


#### **Targeted RNA Therapeutic Modalities**

ModDetect panels are designed to facilitate the detection of specific chemical modifications independent of the sequence or location of the modification and can be used to evaluate a variety of RNA Tx modalities and nucleic acid structures.




#### ModDetect<sup>™</sup> Panels as Analytical Tools

The ModDetect panels serve as analytical tools against different types of nucleic acid targets including modified backbones and 2'-O ribose moieties. Our in-house development ensures supply chain security and respect of client confidentiality and intellectual property.



#### ModDetect<sup>™</sup> PS Panel Results






Figure 2A-B. IF of ModDetect Reagent PS05. Mouse glioma cells derived

from C57 black mice were cultured and treated with Oligo Tx drug. After

fixation with paraformaldehyde, cells were stained with alpha tubulin

at a 1:2000 dilution. Punctate cytoplasmic staining is consistent with

endosomal storage of ASO within the cell, as expected for this Oligo Tx


drug. Vehicle only treated cells showed no staining (not shown).

(red) and ModDetect reagent PS05 (green). Reagent PS05 was used

Figure 1. ELISA of a subset of PS specialty reagents reactive with an FDA-approved gapmer. Streptavidin-coated plates were coated in duplicate with 5 pmol/well of 5'Biotin-MOE-PS Gapmer. The starting dilution was 5 µg/mL and the X-axis represents the Log10 of a 3-fold dilution. This titration is a 4-parameter curve fit where the IC50 is defined as the titer. PS03 IC50: 11 ng/mL [●], PS04 IC50: 9 ng/mL [●], PS05 IC50: 10 ng/mL [▲], PS09 IC50: 18 ng/mL [▼]. Assay performed using Rabbit Anti-Mouse IgG HRP conjugated (#610-403-C46) at 1:8,000 and TMB substrate (#TMBE-1000).

### About the ModDetect<sup>™</sup> Phosphorothioate (PS) Panel

- Includes (5) reagents reactive to phosphorothioate to test and select based on performance with client drug
- Designed for use in various immunoassays, such as ELISA, IF, and IHC
- Includes (3) secondary antibodies
- Unconjugated and biotin-conjugated panels available



Low on time or resources? Our team can perform *in-vitro* assay evaluation studies using our ModDetect panels to determine the best reagent for detecting your oligonucleotide candidates.



**Time Saving** Save 9-12 months in drug development with faster candidate triage



Less Risk Several uniquely-specific reagents are provided for multiple immunoassays



**Cost Saving** Eliminate the need for custom antibody development

