

# BMP Products and Cell Sorting

## I. General Product Description

Biomagnetic Particles (BMP) are a series of solid supports useful in magnetic separation. BMP are ~ 1 um magnetic particles consisting of an iron oxide core with a silane coating. The particle surface is functionalized with amine or carboxyl groups for the covalent attachment of proteins, glycoproteins and other ligands with retention of biological activity. BMP products are also supplied covalently attached to a number of secondary antibodies, such as goat anti-mouse and goat anti-rat as well as streptavidin, biotin and other binding proteins. All BMP solid phases are superparamagnetic, i.e. they respond well to magnetic fields but do not themselves become magnetized. The inability to become permanently magnetized permits repeated magnetic extraction without magnetically induced aggregation. Efficient and rapid removal of BMP particles from suspension is achieved through the application of an external magnetic field.

## II. Principle of Magnetic Cell Sorting

BMP are particularly suitable for cell isolations using either a *direct* or an *indirect* procedure. In the *direct* procedure, monoclonal or polyclonal antibodies to cell surface antigens are covalently attached to BMP before incubation with the cell suspension. In the *indirect* procedure, the cells to be isolated are pretreated with the appropriate monoclonal or polyclonal antibodies and are subsequently magnetically immunoprecipitated with the appropriate magnetic second antibody. Variations of the *indirect* method include pretreating the cells to be isolated with biotin-labeled antibodies and magnetically separating them with BMP Streptavidin. The target cells magnetically labeled by either the *direct* or *indirect* method are easily isolated by applying a rare earth magnet directly against the side of the tube or tissue culture flask such that the cells are magnetically separated perpendicular to gravity. Most cell separations require a 5-10 minute magnetic separation.

## III. Magnetic Cell Sorting Guidelines

No single magnetic cell sorting protocol is applicable for *all* cell separations. The best way to start is by following the general guidelines offered here while referring to published protocols (See References). Each lab should then optimize conditions, i.e. antibody concentrations, incubations times, particle-to-cell ratios to achieve the desired results.

### A. Preparation of Cells

Lymphocytes may be partially purified using a LeucoPREP™ Tube, ficoll gradient, nylon wool, or similar method. However, it may be possible to sort cells directly from whole blood and other sources. Since any particular cell source will have its own unique requirements for purification, procedures must be optimized. Ideally, cell concentrations should **NOT** exceed 5 million total cells per ml. The denser the mixture of cells, the more likely there will be nonspecific binding, clumping and trapping of cells. The investigator should work with the most dilute cellular suspension possible. Sorts of one million cells per ml are most common. If cells clump during washes with sterile media, it may be due to the release of DNA by necrotic cells. These clumps can be easily broken up with the use of 0.1% DNase in the cell medium.

### B. BMP Preparation

BMP products are not claimed to be sterile and are supplied in buffers containing azide. Although not sterile, conjugated BMP products are all prepared by procedures which essentially result in a low bio-burden, i.e. glutaraldehyde conjugations. To remove azide and to prepare the magnetic particles for cell sorting, wash the magnetic particles 2-3 times in an appropriate sterile culture medium (or buffer) containing antibiotics such as penicillin, streptomycin or gentamycin. Washes **MUST** be performed using a BMP magnetic separation unit. **Do not centrifuge during wash steps.**

Media and buffers containing 5-10% protein are recommended to reduce possible non-specific binding. Be advised that too much protein in the medium may inhibit the binding of particles to cells. Proteins such as human serum albumin, bovine serum albumin, fetal bovine serum, or milk solids may be used to lower non-specific binding.

### C. Positive and Negative Selection Using BMP

BMP products may be used in both the *positive* and *negative* selection of cells. Depending upon antigen availability and the size of the target cell population, *positive* and *negative* cell sorting applications may require 20-80 magnetic particles per cell based on the total cell population. Multiple sorts may also be performed for both *positive* and *negative* selection. BMP products contain approximately  $5 \times 10^8$  magnetic particles per mg. Most BMP products are supplied at 1 mg BMP/ml and 5 mg BMP/ml concentrations. Since the particle-to-cell ratio is based upon the **TOTAL** cell population, the following sample calculation applies for both *positive* and *negative* selection.

#### Sample Calculation for Positive or Negative Selection

Consider a system in which there is  $1 \times 10^7$  total cells and in which the target cell population is 30%. In this example we will use a ratio of 50 particles per total cell and the product BMP Goat anti-Mouse IgG (code # 610-119-019) which is supplied as a 1 mg BMP/ml preparation:

$1 \times 10^7$  total cells x 50 particles per total cell =  $5 \times 10^8$  magnetic particles required.

Since 1 mg/ml BMP contain  $5 \times 10^8$  particles per mg, which is equivalent to  $5 \times 10^8$  magnetic particles per ml, the volume of washed BMP required is:

$5 \times 10^8$  magnetic. particles req'd. = 1.0 ml of BMP Goat anti-Mouse IgG

$5 \times 10^8$  magnetic particles per ml

## **D. Incubation Guidelines for Positive and Negative Selection**

1. Incubations should be performed on ice or at 4° C to minimize patching, capping and phagocytosis. Also, cell viability may be best preserved by keeping the cells on ice. However, room temperature or even 37° C may be optimal for certain cell types. The ideal incubation temperature for sorting may vary with the application. Therefore, the investigator may wish to investigate other temperatures should low yield and/or cell viability become concerns.
2. Magnetic particles should be incubated with the cells for 15 to 30 minutes at 4° C. Long incubations are not recommended as magnetic particles may detach from the target cells as a result of cell surface changes over time. During incubation, gentle swirling of the reaction vessel at 10 minute intervals will keep the BMP in suspension. (Continuous rotation is not recommended).
3. To reduce non-specific binding and prevent trapping, cell sorts should be performed in total volumes, which include the cell volume plus the BMP volume, greater than or equal to 1 ml. Additional media or buffer should be added to volumes less than 1 ml to bring the volume to at least 1 ml.

## **E. Separation Guidelines for Positive and Negative Selection**

1. Magnetic separation must be performed perpendicular to gravity with the pellet formed on the side of the flask or tube. This technique is used to keep the unselected cells from contaminating the magnetic pellet due to gravity. The magnet used should have a strength greater than 20 megagauss Orstead. Laboratory Supermagnets™ for use with BMP (code # BMP-00-001 thorough BMP-00-007) have been designed specifically for magnetically separating cells and are recommended.
2. Separation times of 5-10 minutes are generally sufficient for complete separations. A clear supernatant indicates that the separation is complete.
3. Once separation is complete, the supernatant should be removed without disturbing the magnetic pellet. For this reason, vacuum aspiration is not recommended.

## **F. General Procedure for the Indirect Method for Isolating Lymphocytes**

This procedure is for negative or positive cell selection is best carried out in sterile tissue culture tubes, or flasks. The following is an example of a two step protocol which uses BMP Goat anti-Mouse IgG (code # 610-119-019). The target population is assumed to be 20% of the total. All steps are done at 4° C. *The example uses a particle to total cell ratio of 50:1 for BMP Goat anti-Mouse IgG which is supplied at a concentration of 1 mg/ml.*

1. Approximately 10 million total cells are placed in a tube in 10 ml of RPMI with 5% fetal bovine serum and antibiotics.
2. Wash 1 ml of BMP Goat anti-Mouse IgG three times in 1 ml of sterile medium containing antibiotics. Use a magnet to pull the magnetic particles to the side of the tube and shake vigorously to resuspend the magnetic particles during washing. Resuspend in 1 ml of sterile medium.
3. Depending upon the source of the antibody and the manufacturer's recommendations, 5-20 µg of monoclonal antibody per one million target cells is typically needed. Add the appropriate amount of monoclonal antibody to the washed 1 ml of BMP Goat anti-Mouse IgG from step 2 and mix. Incubate at 4° C for 20 minutes.
4. Magnetically separate the BMP Goat anti-Mouse IgG/antibody complex and wash three times with 1 ml of sterile medium. Resuspend in 1 ml of sterile medium.
5. Add the 1 ml of washed BMP Goat anti-Mouse IgG/antibody complex to the 10 ml of cells. (The particle-to-total-cell-ratio is 50:1). Gently swirl the cell/particle mixture to resuspend the cells. The cells are incubated for 20 minutes at 4° C. Swirl the cell/particle suspension every 10 minutes to promote attachment. Magnetically separate for 10 minutes and save the supernatant for a negative selection or save the magnetic pellet for a positive selection.
6. Centrifuge and resuspend the cells in fresh medium for negatively selected cells. For positively selected cells, refer to "Removal of Cells from BMP After Positive Selection, Section H.

### **G. General Procedure for the Direct Method of Separating Lymphocytes**

This procedure is for negative or positive cell selection is best carried out in sterile tissue culture tubes, or flasks. The following is an example of a one step protocol which uses BMP covalently attached to CD8 monoclonal antibody (code # 209-318-354). The target population is assumed to be 25% of the total. All steps are done at 4 ° C. *The example uses a particle to total cell ratio of 100:1 for BMP anti-CD8 which is supplied at a concentration of 1 mg/ml.*

1. Deliver one ml of approximately  $1 \times 10^6$  cells in RPMI with 5% fetal bovine serum and 1% penicillin-streptomycin into an appropriate tube.
2. Wash 0.2 ml of BMP anti-CD8 magnetic particles three times in 0.2 ml of sterile medium containing antibiotics. Use a magnet to pull the magnetic particles to the side of the tube and shake vigorously to resuspend the magnetic particles during washing.
3. Add 0.2 ml of washed BMP anti-CD8 magnetic particles to the cells and swirl to mix the cells with the BMP. (The particle-to-cell-ratio in this example is 100 particles per cell based on the total cell population).
4. Incubate the cells with the BMP anti-CD8 for 20-30 minutes. Swirl the cell/particle suspension once every ten minutes to promote attachment. Magnetically separate the supernatant for 10 minutes twice and save the supernatant for a negative selection or save the magnetic pellet for a positive selection.
5. Centrifuge and resuspend the negatively selected cells in fresh medium for negative selection. For positively selected cells, refer to "Removal of Cells from BMP after Positive Selection, Section H.

#### **H. Removal of Cells from BMP after Positive Selection**

1. Methods for detaching magnetic particles from cells after separation include culturing cells for up to 48 hours during which magnetic particles fall away from the cells due to cell surface turnover (See References, Pricop et. al.) or using a protease, such as chymopapain, to break the antigen-antibody bond. Each of these procedures has limitations. Not all particles may detach from cells during culturing and the use of a protease may damage cells. Depending upon the application, it may **NOT** be necessary to remove the cells from the BMP. BMP are only 1  $\mu\text{m}$  in size and can be successfully used in flow cytometry equipment as they will not jam the equipment and are distinguishable from cells. Alternatively, negative selection should be considered.

#### **IV. Magnetic Cell Sorting References**

Adkins, B., Ghanei, A., and Hamilton, K. 1993. Developmental Regulation of IL-4, IL-2, and IFN- $\delta$  Production by Murine Peripheral T Lymphocytes. Journal of Immunology 151:6617-6626.

Adkins, B. and Hamilton, K. 1992. Freshly Isolated, Murine Neonatal T Cells Produce IL-4 in Response to anti-CD3 Stimulation. Journal of Immunology 149:3448-3455.

Agematsu, K., Kobata, T., Sugita, K., Freeman, G., Beckmann M., Schlossman, S. and Morimoto, C. 1994. Role of CD27 in T Cell Immune Response. Journal of Immunology 153:1421-1429.

Agematsu, K., Kobata, T., Sugita, K., Hirose, T., Schlossman, S., and Morimoto, C. 1995. Direct Cellular Communications Between CD45RO and CD45RA T Cell Subsets Via CD27/CD70. Journal of Immunology 154:3627-3635.

Alzona, M., Jack, H.M., Fisher, R.I. and Ellis, T.M. 1995. IL-12 Activates IFN- $\gamma$  Production Through the Preferential Activation of CD30+ T Cells. Journal of Immunology 154:9-16.

Askenase, P., Szczepanik, M., Ptak, M., Paliwal, V., and Wlodzimierz, P. 1995.  $\gamma$   $\delta$  T Cells in Normal Spleen Assist Immunized  $\alpha$   $\beta$  T Cells in the Adoptive Cell Transfer of Contact Sensitivity. Journal of Immunology 154:3644-3653.

Bathon, Joan M., Croughan, John C., MacGlashan, Donald W., and Proud, David. 1994. Bradykinin Is A Potent and Relatively Selective Stimulus for Cytosolic Calcium Elevation in Human Synovial Cells. Journal of Immunology 153:2600-2608.

Beno, D., Stover, A., and Mathews, H. 1995. Growth Inhibition of *Candida albicans* Hyphae by CD8+ Lymphocytes. Journal of Immunology 154:5273-5281.

Berman, Jeffrey S., Center, David M. 1987. Chemotactic Activity of Porcine Insulin for Human T Lymphocytes in Vitro. Journal of Immunology 138:2100.

Bieva, Claude J., VanderBruggen, F.J., Stryckmans, P.A. 1989. Malignant Leukemic Cell Separation by Iron Colloidal Immunomagnetic Adsorption. Experimental Hematology 17:914-920.

Bonomo, A., Kehn, P., Payer, E., Rizzo, L., Cheever, A., and Shevach, E. 1995. Pathogenesis of Post-Thymectomy Autoimmunity : Role of Syngeneic MLR-Reactive T Cells. Journal of Immunology 154:6602-6611.

Brissette-Storkus, C. et al. 1994. Characterization and Function of the NKR-P1dim/T Cell Receptor- $\alpha$   $\beta$  + Subset of Rat T Cells. Journal of Immunology 152:388-396.

Chen, Y., Kuchroo, V.K., Inobe, J., Hafler, D. and Weiner, H.L. 1994. Regulatory T Cell Clones Induced by Oral Tolerance: Suppression of Autoimmune Encephalomyelitis. Science 265:1237-1240.

Chervenak, R., Dempsey, D., Soloff, R., Wolcott, M. and Jennings, S. 1993. The Expression of CD4 by T Cell Precursors Resident in Both the Thymus and the Bone Marrow. Journal of Immunology 151:4486-4493.

Chervonsky, Alexander V., Golovkina, Tatyana V., Ross, Susan R., and Janeway, Charles A. 1995. Differences in the Avidity of TCR Interactions with a Superantigenic Ligand Affect Negative Selection

but do not Allow Positive Selection. Journal of Immunology 155:5115-5123.

Chilton, P. and Fernandez-Botran, R. 1993. Production of Soluble IL-4 Receptors by Murine Spleen Cells is Regulated by T Cell Activation and IL-4. Journal of Immunology 151:5907-5917.

Conlon, K., Ochoa, A., Kopp, W., Ortaldo, J., Urba, W., Longo, D., and Young, H. 1992. Enhanced Lymphokine Production and Lymphokine Receptor Expression in Multiple Antibody-Stimulated Human CD4+ Peripheral Blood Lymphocytes. Journal of Immunology 149:3278-3289.

Cruikshank, William W., Berman, Jeffrey S., Theodore, Arthur C., Bernardo, John, and Center, David M. 1987. Lymphokine Activation of T4 and T Lymphocytes and Monocytes. Journal of Immunology 138:3817-3823.

Dieli, F., Asherson, G., Romano, G., Sireci, G., Gervasi, F., and Salerno, A. 1994. IL-4 is Essential for the Systemic Transfer of Delayed Hypersensitivity by T Cell Lines. Journal of Immunology 152:2698-2704.

Dong, Z. and Wortis, H. 1994. Function of Bone Marrow Stromal Cell Lines Derived from Nude Mice. Journal of Immunology 153:1441-1448.

Elson, C.O., Holland, S.P., Dertzbaugh, M.T., Cuff, C.R. and Anderson, A.O. 1995. Morphologic and Functional Alterations of Mucosal T Cells by Cholera Toxin and its B Subunit. Journal of Immunology 154:1032-1040.

Ettinger R. et al. 1994. Functional Distinctions between MRL-*lpr* and MRL-*gld* Lymphocytes. Journal of Immunology 152:1557-1568.

Ferbas, John J., Toso, John R., Logar, Alison J., Navratil, Jeannine S., and Rinaldo, Charles R. 1994. CD4+ Blood Dendritic Cells Are Potent Producers of IFN- $\alpha$  in Response to In Vitro HIV-1 Infection. Journal of Immunology 152:4649-4662.

Ford, John R. and Terzaghi-Howe, M. 1992. Characteristics of Magnetically Separated Rat Tracheal Epithelial Cell Populations. Am. J. Physiol. 263:L568-L574.

Foroozan, M., Shi, Y., Shirazian, D., Morgante, L., Miller, A. and Grob, D. 1994. Defective Production of Anti-Inflammatory Cytokine, TGF- $\beta$ , by T Cell Lines of Patients with Active Multiple Sclerosis. Journal of Immunology 152:6003-6010.

Fowler, D., Kurasawa, K., Husebekk, A., Cohen, P. and Gress, R. 1994. Cells of Th2 Cytokine Phenotype Prevent LPS-Induced Lethality during Murine Graft-versus-Host Reaction. Journal of Immunology 152:1004-1013.

Freedman, Mark S., Blain, Manon, and Antel, Jack P. 1991. Differential Responses of CD4+CD45RA+ and CD4+CD29+ Subsets to Activated CD8+ Cells: Enhanced Stimulation of the CD4+CD45RA+

Subset by Cells from Patients with Multiple Sclerosis. Cellular Immunology 133:306-316.

Gilat, D., Hershkovich, R., Mekori, Y. A., Vlodaysky, I. and Lider O. 1994. Regulation of Adhesion of CD4+ T Lymphocytes to Intact or Heparinase-Treated Subendothelial Extracellular Matrix by Diffusible or Anchored RANTES and MIP-1 $\beta$ . Journal of Immunology 153:4899-4906.

Gonzalez, M., Merino, R., Gonzalez, A., and Merino, J. 1995. The Ability of B Cells to Participate in Allogeneic Cognate T-B Cell Interactions In Vitro Depends on the Presence of CD4+ T Cells During Their Development. Journal of Immunology 155:1091-1100.

Gribben, John G. et al. 1992. Bone Marrows of Non-Hodgkin's Lymphoma Patients With a bcl-2 Translocation Can Be Purged of Polymerase Chain Reaction-Detectable Lymphoma Cells Using Monoclonal Antibodies and Immunomagnetic Bead Depletion. Blood 80:1083-1089.

Gryllis, Chryssa, Wainberg, Mark A., Gornitsky, Mervyn and Brenner, Bluma. 1990. Diminution of Inducible Lymphokine-Activated Killer Cell Activity in Individuals with AIDS-Related Disorders. AIDS 4:1205-1212.

Haregewoin, Abebe, Soman, Gopalan, Hom, Richard C., Finberg, Robert W. 1989. Human T Cells Respond to Mycobacterial Heat Shock Protein. Nature 340:309-312.

Hermentin, P. et al. 1990. Hinge-Thiol Coupling of Monoclonal Antibody to Silanized Iron Oxide Particles and Evaluation of Magnetic Cell Depletion. Bioconjugate Chemistry 1:411-418.

Johnston, J., Taub, D., Lloyd, A., Conlon, K., Oppenheim, J. and Kevlin, D. 1994. Human T Lymphocyte Chemotaxis and Adhesion Induced by Vasoactive Intestinal Peptide. Journal of Immunology 153:1762-1768.

June, Carl H., Ledbetter, Jeffrey A., Gillespie, Marjorie M., Lindsten, Tullia and Thompson, Craig B. 1987. T Cell Proliferation Involving the CD28 Pathway is Associated with Cyclosporine-Resistant Interleukin-2 Gene Expression. Molecular and Cell Biology 7:4472-4480.

Kobata, T., Agematsu, K., Kameoka, J., Schlossman, S., and Morimoto, C. 1994. CD27 is a Signal-Transducing Molecule Involved in CD45RA+ Naive T Cell Costimulation. Journal of Immunology 153:5422-5432

Kornfeld, Hardy, Cruikshank, William W., Pyle, Stephen W., Berman, Jeffrey S., Center, David M. 1988. Lymphocyte Activation by HIV-1 Envelope Glycoprotein. Nature 335:445-448.

Laberge, S., Cruikshank, W.W., Beer, D.J., and Center, D.M. 1996. Secretion of IL-16 (Lymphocyte Chemoattractant Factor) from Serotonin-Stimulated CD8+ Cells In Vitro. Journal of Immunology 156:310-315.

Lagoo, A. et al. 1994. IL-2, IL-4, and IFN- $\gamma$  Gene Expression Versus Secretion in Superantigen-

Activated T Cells. Journal of Immunology 152:1641-1652.

Lederman S. et al. 1994. T-BAM/CD40-L on Helper T Lymphocytes Augments Lymphokine-Induced B Cell Ig Isotype Switch Recombination and Rescues B Cells from Programmed Cell Death. Journal of Immunology 152:2163-2171.

Lee, William T. and Vitetta, Ellen S. 1991. The Differential Expression of Homing and Adhesion Molecules on Virgin and Memory T Cells in the Mouse. Cellular Immunology 132:215-222.

Levin, D. Constant, S., Pasqualini, T., Flavell, R., and Bottomly, K. 1993. Role of Dendritic Cells in the Priming of CD4 T Lymphocytes to Peptide Antigen In Vivo. Journal of Immunology 151:6742-6750.

Linette, Gerald P., Hartzman, Robert J., Ledbetter, Jeffrey A., June, Carl H. 1988. HIV-1-Infected Cells Show A Selective Signaling Defect After Perturbation of CD3-Antigen Receptor. Science 241:573.

Loeffler, Cynthia M. et al. 1991. Antitumor Effects of Interleukin 2 Liposomes and Anti-CD3-Stimulated T Cells Against Murine MCA-38 Hepatic Metastasis. Cancer Research 51:2127-2132.

Lopes, M.G., daVeiga, V.F., Santos, A.R., Fonseca, M.E., DosReis, G.A. 1995. Activation-Induced CD4+ T Cell Death by Apoptosis in Experimental Chagas' Disease. Journal of Immunology 154:744-752.

Mehrotra, P., Grant, A., and Siegel, J. 1995. Synergistic Effects of IL-7 and IL-12 on Human T Cell Activation. Journal of Immunology 154:5093-5102.

Moran, T., Toy, E., Kuzu, Y., Kuzu, H., Isobe, H., and Schulman, J. 1995. Staphylococcal Enterotoxin B-Activated T Cells Can Be Redirected to Inhibit Multicycle Virus Replication. Journal of Immunology 155:759-765.

Nagata, M. Santamaria, P. Kawamura, T. Utsugi, T. and Yoon, J. 1994. Evidence for the Role of CD8 + Cytotoxic T Cells in the Destruction of Pancreatic  $\beta$  -Cells in Nonobese Diabetic Mice. 1994. Journal of Immunology 152:2042-2050.

Nestle, F., Zheng, X., Thompson, C., Turka, L., and Nickoloff, B. 1993. Characterization of Dermal Dendritic Cells Obtained from Normal Human Skin Reveals Phenotypic and Functionally Distinctive Subsets. Journal of Immunology 151:6535-6545.

Padmanabhan, Raji, Corsico, Christopher D., Howard, Tazuko H., Holter, W., Fordis, C., Willingham, Mark, and Howard, Bruce. 1989. Purification of Transiently Transfected Cells by Magnetic Affinity Cell Sorting. Analytical Biochemistry 170:341-348.

Padmanabhan R., Corsico, C., Holter, W., Howard, T. and Howard B.H. 1989. Purification of Transiently Transfected Cells by Magnetic-Affinity Cell Sorting. Journal of Immunogenetics 116:91-

102.

- Pilaro, Anne M., Taub, Dennis D., McCormick, Kathryn L., Williams, Hope M., Sayers, Thomas J., Fogler, William E., Wiltrout, Robert H. 1994. TNF- $\alpha$  Is a Principal Cytokine Involved in the Recruitment of NK Cells to Liver Parenchyma. Journal of Immunology **153**:333-342.
- Pricop, L. et al. 1993. Characterization of the Fc $\mu$  Receptor on Human Natural Killer Cells. Journal of Immunology **151**:3018-3029.
- Rabinowich H., Lin, W., Herberman, R., and Whiteside T. 1994. Signaling Via CD7 Molecules on Human NK Cells. Journal of Immunology **153**:3504-3513.
- Rabinowich H., Sedlmayr, P., Herberman, R., and Whiteside, T. 1993. Response of Human NK Cells to IL-6 Alterations of the Cell Surface Phenotype, Adhesion to Fibronectin and Laminin, and Tumor Necrosis Factor- $\alpha$  / $\beta$  Secretion. Journal of Immunology **150**:4844-4855.
- Rabinowich, Hannah et al. 1992. Role of Cytokines in the Adoptive Immunotherapy of an Experimental Model of Human Head and Neck Cancer by Human IL-2 Activated Natural Killer Cells. Journal of Immunology **149**:340-349.
- Rabinowich, H., Lin, W., Amoscato, A., Herberman, R.B. and Whiteside, T.L. 1995. Expression of Vitronectin Receptor on Human NK Cells and Its Role in Protein Phosphorylation, Cytokine Production and Cell Proliferation. Journal of Immunology **154**:1124-1135.
- Rasmussen, R., Counts, S., Daley, J. and Schlossman, S. 1994. Isolation and Characterization of CD6- T Cells from Peripheral Blood. Journal of Immunology **152**:527-536.
- Rosenman, S., Ganji, A., Tedder, T., and Gallatin, W. 1993. Syn-capping of human T Lymphocyte Adhesion/Activation Molecules and their Redistribution During Interaction with Endothelial Cells. Journal of Leukocyte Biology **53**:1-10.
- Saha, K., and Wong, P. 1992. Ts1, a Temperature-Sensitive Mutant of Moloney Murine Leukemia Virus TB, Can Infect Both CD4+ and CD8+ T Cells but Requires CD4+ T Cells in Order to Cause Paralysis and Immunodeficiency. Journal of Virology **66**:2639-2646.
- Saitoh, S-I., Kosugi, A., Noda. S., Yamamoto, N., Ogata, O., Minami, Y., Miyake, K., and Hamaoka, T. 1995. Modulation of TCR-Mediated Signaling Pathway by Thymic Shared Antigen-1 (TSA-1)/Stem Cell Antigen-2 (Sca-2). Journal of Immunology **155**:5574-5581.
- Santamaria, P., Utsugi, T., Park, B., Averill, N., Kawazu, S., and Yoon, J. 1995. Beta-Cell-Cytotoxic CD8+ T Cells from Nonobese Diabetic Mice Use Highly Homologous T. Cell Recetor  $\alpha$  -Chain CDR3 Sequences. Journal of Immunology **154**:2494-2503.
- Sen-Majudmar, A., Guidos, C., Kina, T. Lieberman, M. and Weissman, I. 1994. Characterization of

Preneoplastic Thymocytes and of Their Neoplastic Progression in Irradiated C57BL/Ka Mice. Journal of Immunology 153:1581-1592.

Selvey, L.A., Morse, H.C., June, C.H. and Hodes, R.J. 1995. Analysis of Antigen Receptor Signaling in B Cells from Mice with a Retrovirus-Induced Acquired Immunodeficiency Syndrome. Journal of Immunology 154:171-179.

Shen, F. Xu, X.L., Graf, L.H. and Chong, A.S.-F. 1995. CD45-Cross-Linking Stimulates IFN- $\gamma$  Production in NK Cells. Journal of Immunology 154:644-652.

Smyth, Mark J., Strobl, Susan L., Young, Howard A., Ortaldo, John R. and Ochoa, Augusto C. 1991. Regulation of Lymphokine-Activated Killer Activity and Pore Forming Protein Gene Expression in Human Peripheral Blood CD8+ T Lymphocytes. Journal of Immunology 146:3289-3297.

Spangrude, Gerald J and Scollay, Roland. 1990. A Simplified Method for Enrichment of Mouse Hematopoietic Stem Cells. Experimental Hematology 18:920-926.

Spangrude, Gerald J., Heimfeld, Shelly and Weissman, Irving L. 1988. Purification and Characterization of Mouse Hematopoietic Stem Cells. Science 241:58-62.

Stout, R.D., Suttles, J., Xu, J., Grewal, I.S. and Flavell, R.A. 1996. Impaired T Cell-Mediated Macrophage Activation in CD40 Ligand-Deficient Mice. Journal of Immunology 156:8-11.

Sugie, K., Minami, Y., Kawakami, T., and Uchida, A., 1995. Stimulation of NK-Like YT Cells via Leukocyte Function-Associated Antigen (LFA)-1. Journal of Immunology 154:1691-1698.

Tacke, M., Clark, G., Dallman, M., and Hunig T. 1995. Cellular Distribution and Costimulatory Function of Rat CD28. Journal of Immunology 154:5121-5127.

Taylor, R., Pocanic, F., Reist, C., and Wright, E. 1991. Complement-Opsonized IgG Antibody/dsDNA Immune Complexes Bind to CR1 Clusters on Isolated Human Erythrocytes. Clinical Immunology and Immunopathology 61:143-160.

Trickett, A.E., Ford, D.J., Lam-Po-Tang, P.R.L., and Vowels, M.R. 1990. Comparison of Magnetic Particles for Immunomagnetic Bone Marrow Purging Using an Acute Lymphoblastic Leukaemia Model. Transplantation Proceedings 22:2177-2178.

Turka, L., Schatz, D., Oettinger, M., Chun, J., Gorka, C., Lee, K., McCormack, W., and Thompson, C. 1991. Thymocyte Expression of RAG-1 and RAG-2: Termination by T Cell Receptor Cross-Linking. Science 253:778-781.

Vaccaro, Dennis E. 1990. Applications of Magnetic Separation: Cell Sorting. American Biotechnology Laboratory. April 1990.

Vujanovic, N.L., Yashumura, S., Hirabayashi, H., Lin, W.C., Watkins, S. Herberman, R.B. and Whiteside, T.L. 1995. Antitumor Activities of Subsets of Human IL-2-Activated Natural Killer Cells in Solid Tissues. Journal of Immunology 154:281-289.

Wirth, S., Guidotti, L., Ando, K., Schlicht, H. and Chisari, F. 1995. Breaking Tolerance Leads to Autoantibody Production but Not Autoimmune Disease in Hepatitis B Virus Envelope Transgenic Mice. Journal of Immunology 154:2504-2515.

Witkowski, Jacek M., Li, Shaokang P., Gorgas, G., and Miller, Richard A. 1994. Extrusion of the P Glycoprotein Substrate Rhodamine-123 Distinguishes CD4 Memory T Cell Subsets That Differ in IL-2 Driven IL-4 Production. Journal of Immunology 153:658-665.

Yabkowitz, R., Dixit, V., Guo, N., Roberts, D. and Shimizu, Y. 1993. Activated T-Cell Adhesion to Thrombospondin is Mediated by the  $\alpha 4\beta 1$  (VLA-4) and  $\alpha 5\beta 1$  (VLA-5) Integrins. Journal of Immunology 151:149-158.

Yuen, M-H., Protti, M., Diethelm-Okita, B., Moiola, L., Howard, J., Conti-Fine, B. 1995. Immunoregulatory CD8+ Cells Recognize Antigen-Activated CD4+ Cells in Myasthenia Gravis Patients and in Healthy Controls. Journal of Immunology 154:1508-1520.

Additional information on BMP and BMP antibody conjugates is available from our Technical Service staff upon request. Please call 8 a.m. to 5 p.m. EST, fax or e-mail your request at any time. Access this document from the web at [www.rockland-inc.com/morebmpinfo.html](http://www.rockland-inc.com/morebmpinfo.html).