

Product Specification Sheet

Product: anti-Chicken Red Blood Cells (RBC) [Rabbit]

Code: 103-4139

Lot # 1204

Size: 2.0 ml

Physical State: Lyophilized

Protein Concentration: 80.0 mg/ml (by Refractometry)

Buffer: 0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.4

Stabilizer: None

Preservative: None

Application(s): suitable for agglutination and related assays.

Recommended Dilution(s): An agglutination titer against serial dilutions of washed pooled cells on a titer plate of 1:64 was observed for this product. Optimal titers for applications should be determined by the researcher.

Storage Conditions: Store vial at 4° C prior to restoration. Restore with 2.0 ml of deionized water (or equivalent). Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.

Purity: This product was prepared from polyspecific antiserum by a delipidation and defibrination.

Immunogen: Chicken washed pooled Red Blood Cells (RBC)

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.