

Product Specification Sheet

Product: Mouse Liver powder

Code: C118

Lot # 4008

Size: 1.0 g

Physical State: Powder

Buffer: None

Stabilizer: None

Preservative: None

Preparation Method: Tissue powders are produced by washing the various organs repeatedly in saline to remove blood followed by multiple precipitations in acetone until all lipid is removed. Acetone is removed by desiccation.

Recommended Use: Suitable for use as a specific adsorbent to remove unwanted reactivities or as a source material. Product provides a high nuclear content and is micronized to a fine powder.

Storage Conditions: Store vial at 2 to 8° C prior to opening. No expiration date is given for this product if properly stored.

Strain: Swiss-Webster

Sex: Mixed

USDA Certification: All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation.

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.