

Certificate of Analysis

Product: Affinity Purified Anti-apoLipoprotein C-I [Goat] Minimum Cross Reactivity to Type A-I, A-II, B, C-II, C-III, E and J apoLipoproteins.

Code: 600-101-112

Lot # 18864

Size: 1.0 mg

Physical State: Liquid (sterile filtered)

Antibody Concentration: 1.0 mg/ml (by UV absorbance at 280 nm)

Buffer: 0.125M Sodium Borate, 0.075M Sodium Chloride, 0.005M EDTA; pH 8.0

Stabilizer: None

Preservative: 0.01% (w/v) Sodium Azide

Application(s): Anti-apoLipoprotein antibodies have been used for indirect trapping ELISA for quantitation of antigen in serum using a standard curve, for immunoprecipitation and for western blotting for highly sensitive qualitative analysis. This product has been tested by immunohistochemistry using paraffin embedded tissues.

Recommended Dilution(s):	ELISA	1:5,000 - 1:10,000
	WESTERN BLOT	1:5,000 - 1:10,000
	IMMUNOHISTOCHEMISTRY	1:50 - 1:200
	OTHER APPLICATIONS	User Optimized

Storage Conditions: Store vial at 4° C prior to opening. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, mix with an equal volume of glycerol, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of opening product.

Purity and Specificity: This product has been prepared by immunoaffinity chromatography using immobilized antigens followed by extensive cross-adsorption against other apoLipoproteins and human serum proteins to remove any unwanted specificities. Typically less than 1% cross reactivity against other types of apoLipoprotein was detected by ELISA against purified standards. This antibody reacts with human apoLipoprotein C-I and has negligible cross-reactivity with Type A-I, A-II, B, C-II, C-III, E and J apoLipoproteins. Specific cross reaction of anti-apoLipoprotein antibodies with antigens from other species has not been determined. Non-specific cross reaction of anti-apoLipoprotein antibodies with other human serum proteins is negligible.

Immunogen: apoLipoprotein Type C-I produced synthetically in full-length form (not selected epitopes) using conventional peptide technology.

USDA Certification: All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation.

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.