

Product Specification Sheet

Product: Guinea Pig Complement (fresh frozen)

Code: C300-0010

Lot # 19093

Size(s): 10 ml

Physical State: Fresh Frozen

Protein Concentration: 79.0 mg/ml (by Refractometry)

IgG Concentration: Not determined

Buffer: None

Stabilizer: None

Preservative: None

Recommended Dilution(s): Not specified.

Storage Conditions: Store vial at -70° C prior to opening. **NOTE: COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY.** Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of opening.

Immunoelectrophoresis: Normal

pH: Normal

Hemoglobin: Normal

Strain: Mixed breed (males only)

Sex: Male

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.