



800-656-7625

fax 610-367-7825

Product Specification Sheet

Product: *MaxTag™ Histo* 30% (w/v) Hydrogen Peroxide Solution

Code: KHJ001

Lot # 16283

Size: 10 ml

Physical State: Clear colorless liquid

Molecular Formula: H₂O₂

Molecular Weight: 34.01

Stabilizer: None

Preservative: None

Application(s): Please read the instruction insert for the *MaxTag™ Histo* kit.

Storage Conditions: 2° to 8° C. Keep dry. Warm to room temperature before opening.

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.