



800-656-7625

fax 610-367-7825

Certificate of Analysis

Product: Affinity Purified Anti-Human FANCE [Rabbit]

Code: 600-401-673

Lot #: 21603

Size: 100 µg

Physical State: Liquid (sterile filtered)

Antibody Concentration: 1.0 mg/ml (by UV absorbance at 280 nm)

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Stabilizer: None

Preservative: 0.01% (w/v) Sodium Azide

Storage Conditions: Store vial at -20° C prior to opening. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of opening.

Application Note(s): This affinity purified antibody has been tested for use in ELISA against the immunizing peptide. Reactivity in other immunoassays is unknown.

Recommended Dilutions: This product has been assayed by ELISA against 0.1 µg of the immunizing peptide. A 1:4,000 to 1:18,000 dilution of the antibody is recommended for this assay.

Purity and Specificity: This is an affinity purified antibody produced by immunoaffinity chromatography using the immunizing peptide after immobilization to a solid phase.

Immunogen: This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to amino acids 1-13 of Human FANCE. This Fanconi anemia protein promotes the nuclear accumulation of FANCC. FANCE protein is part of the FA nuclear complex, binding FANCC & FANCD2. It is required for the nuclear accumulation of FANCC and provides a bridge between the FA complex and FANCD2. FANCC mutants do not bind FANCE or prevent chromosome breakage. LocusLink ID 2178; SwissProt accession number Q9HB96

Peptide sequence: M-A-T-P-D-A-G-L-P-G-A-E-G-C

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 326, Gilbertsville, Pennsylvania, USA.