

Certificate of Analysis

Product: DyLight™ 800 Conjugated Affinity Purified anti-Mouse IgG (H&L) [Goat] Minimum Cross Reactivity to Bovine, Chicken, Goat, Guinea Pig, Hamster, Horse, Human, Rabbit, Rat and Sheep Serum Proteins

Code: 610-145-121

Lot # 20314

Size: 100 µg

Physical State: Lyophilized

Antibody Concentration: 1.0 mg/ml (by UV absorbance at 280 nm)

Label: DyLight™ 800 (MW 1050)

Fluorochrome/Protein Ratio: 2.0 DyLight™ 800 per mole of Goat IgG

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Excitation Wavelength:	Emission Wavelength:
770 nm (in PBS)	794 nm (in PBS)

Stabilizer: 10 mg/ml BSA, IgG and Protease free

Preservative: 0.01% (w/v) Sodium Azide

Application Note(s): This antibody is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This antibody is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

Instrument compatibility: The emission spectra for this DyLight™ conjugate match the principal output wavelengths of most common fluorescence instrumentation.

Recommended Dilution(s):	FLUORESCENT WESTERN BLOT	>1:20,000
	FLISA	>1:10,000
	IMMUNOFLUORESCENT MICROSCOPY	>1:5,000
	OTHER APPLICATIONS	User Optimized

Storage Conditions: Store vial at 4° C prior to reconstitution. Rehydrate with 0.1 ml of deionized water (or equivalent). This product is stable for several weeks at 4° C as an undiluted liquid. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use. Expiration date is one (1) year from date of receipt if properly stored.

Purity and Specificity: This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Mouse IgG and Mouse Serum. No reaction was observed against Bovine, Chicken, Goat, Guinea Pig, Hamster, Horse, Human, Rabbit, Rat and Sheep Serum Proteins. This antibody will react with heavy chains of mouse IgG and with light chains of most mouse immunoglobulins.


Immunogen: Mouse IgG, whole molecule






Conjugation Chemistry: N-hydroxysuccinimide (NHS) ester

Properties of DyLight™ Fluorescent Dyes:

e-mail: info@rockland-inc.com

Visit our website at: www.rockland-inc.com




Emission	Color	DyLight™ Dye	Ex/Em (nm)	ϵ (M ⁻¹ cm ⁻¹)	Similar Dyes
Green		488	493/518	70,000	Alexa™ 488, Cy2®, FITC
Yellow		549	550/568	150,000	Alexa™ 546, Alexa 555, Cy3®, TRITC
Red		649	646/674	250,000	Alexa™ 647, Cy5®
Near Infrared		680	682/715	140,000	Alexa™ 680, Cy5.5®, IRDye™ 700
Infrared		800	770/794	270,000	IRDye™ 800

Relevant Links: See www.rockland-inc.com/DyLight for additional information. Assortments of DyLight™ conjugates can be ordered as 3 and 10 vial value packs.

Related Products:

- #610-141-121 DyLight™ 488 Conjugated Affinity Purified anti-MOUSE IgG (H&L) [GOAT] MX10
- #610-142-121 DyLight™ 549 Conjugated Affinity Purified anti-MOUSE IgG (H&L) [GOAT] MX10
- #610-143-121 DyLight™ 649 Conjugated Affinity Purified anti-MOUSE IgG (H&L) [GOAT] MX10
- #610-144-121 DyLight™ 680 Conjugated Affinity Purified anti-MOUSE IgG (H&L) [GOAT] MX10
- #610-145-121 DyLight™ 800 Conjugated Affinity Purified anti-MOUSE IgG (H&L) [GOAT] MX10
- #MB-070 Blocking Buffer for Fluorescent Western Blotting

DyLight™ 800 Fluorescence Spectra:



Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. There is no implied license hereunder for any commercial use. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 326, Gilbertsville, Pennsylvania, USA.

DyLight™ is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. IRDye™ and Odyssey® are trademarks of LI-COR Inc. Cy2®, Cy3®, Cy5® and Cy5.5® are trademarks of GE Healthcare Inc. Alexa® is a trademark of Invitrogen Inc.