

Product Specification Sheet

Product: Agarose Immobilized IgG Fraction of anti-GST [Goat]

Code: 200-119-200

Lot # 54123

Size: 20.0 mg antibody

Physical State: Suspension of agarose beads

Antibody/Agarose: 10.0 mg antibody per cc agarose

Agarose/Suspension: 10.0 cc agarose per ml of suspension

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Stabilizer: None

Preservative: 0.01% (w/v) Sodium Azide

Application(s): Suitable for molecular biology, enzyme purification or to remove unwanted cross-reactivity.

Recommended Dilution(s): Not specified. Optimal titers for applications should be determined by the researcher.

Storage Conditions: Store vial at 4° C prior to opening. DO NOT FREEZE. Expiration date is one (1) year from date of opening.

Purity: This product is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of GST or GST containing fusion protein.

Immunogen: Glutathione-S-Transferase (*Schistosoma japonicum*)

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.