

**Product Specification Sheet**

**Product:** Non-Sterile Bovine/Calf Plasma in EDTA, Sodium

**Code:** D400-04-0050

**Lot #** 14921

**Size(s):** 50 ml

**Physical State:** Liquid

**Total Protein Concentration:** 86.0 mg/ml (by Refractometry)

**IgG Concentration:** Normal

**Stabilizer:** None

**Preservative:** None

**Recommended Dilution(s):** Not specified.

**Storage Conditions:** Store vial at -20° C prior to opening. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of opening.

**Sterilization:** Not performed.

**Immunoelectrophoresis:** Normal

**pH:** Normal

**Strain:** Mixed breed and sex

**USDA Certification:** All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation.

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.