

Product Specification Sheet

Product: *MaxTag™ Histo* DAB (3,3'-Diaminobenzidine) Substrate Tablets

Code: DAB-525

Lot: 14259

Size : 5 X 2.5 mg

Physical State: tablet

Applications: This product is designed to be used with the *MaxTag™ Histo* Immunohistochemistry kit. Dissolve one tablet in 5 ml of deionized water. Use liquid substrate within one hour. Immediately before use add hydrogen peroxide (H₂O₂) to 0.1% (v/v).

Storage Conditions: Store container at -20° C prior to opening. Keep dry. Warm to room temperature before opening. Expiration date is one (1) year from date of receipt.

Shipping Conditions: Use care. Ship as a hazardous material. DAB is a suspected carcinogen. Avoid exposure to eyes, skin and mucous membranes. Do not inhale or ingest. Wash hands thoroughly after using.

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.