

**Certificate of Analysis****Product:** Anti-Biotin [Goat]**Code:** 100-1198**Lot #** 2170**Size:** 2.0 ml**Physical State:** Lyophilized**Protein Concentration:** 85.0 mg/ml (by Refractometry)**Buffer:** 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2**Stabilizer:** None**Preservative:** 0.01% (w/v) Sodium Azide**Application(s):** Suitable for immunoblotting (western or dot blot), ELISA, immunoprecipitation and most immunological methods requiring high titer and specificity.

<b>Recommended Dilution(s):</b>	<b>ELISA</b>	1:100,000 - 1:500,000
	<b>WESTERN BLOT</b>	1:10,000 - 1:50,000
	<b>IMMUNOHISTOCHEMISTRY</b>	User Optimized
	<b>OTHER APPLICATIONS</b>	User Optimized

**Storage Conditions:** Store vial at 4° C prior to restoration. Restore with 2.0 ml of deionized water (or equivalent). Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.**Purity and Specificity:** This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against biotinylated IgG and Biotinylated Bovine Serum Albumin (BSA).**Immunogen:** Biotin conjugated to Keyhole Limpet Hemocyanin (KLH)**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 326, Gilbertsville, Pennsylvania, USA.