

Certificate of Analysis

Product: *MaxTag*[™] DAB (3,3'-Diaminobenzidine) Substrate Tablets

Code: DAB-50

Lot : 19904

Size : 50 tablets (5 mg each)

Physical State: Tablet

Applications: This product is designed to be used with the *MaxTag*[™] Immunoblotting kit. Dissolve one tablet in 10 ml of deionized water. Use liquid substrate within one hour. Immediately before use add hydrogen peroxide (H₂O₂) to 0.1% (v/v).

Storage Conditions: Store container at -20°C prior to opening. Keep dry. Warm to room temperature before opening. Expiration date is one (1) year from date of receipt.

Shipping Conditions: Use care. Ship as a hazardous material. DAB is a suspected carcinogen. Avoid exposure to eyes, skin and mucous membranes. Do not inhale or ingest. Wash hands thoroughly after using.

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.

All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation.

All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 326, Gilbertsville, Pennsylvania, USA.