

**Product Specification Sheet**

**Product:** Guinea Pig Complement (lyophilized)

**Code:** C200-0005

**Lot #** 12537

**Size:** 5.0 ml

**Physical State:** Lyophilized

**Protein Concentration:** 85.0 mg/ml (by Refractometry)

**IgG Concentration:** Not determined

**Buffer:** None

**Stabilizer:** None

**Preservative:** None

**Recommended Dilution(s):** User optimized.

**Storage Conditions:** Store vial at 4° C prior to restoration. Restore with 5.0 ml of diluent (supplied with product). Centrifuge product (at 2° to 8° C) if not completely clear. Aliquot contents and freeze at -70° C or below. **NOTE: COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY.** Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.

**Immunoelectrophoresis:** Normal

**pH:** Normal

**Hemoglobin:** Normal

**Strain:** Mixed breed

**Sex:** Male

**USDA Certification:** All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation.

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.