

**Certificate of Analysis**

**Product:** Rat Gamma Globulin

**Code:** D611-0050

**Lot #** 19617

**Size:** 50.0 mg

**Physical State:** Lyophilized

**Antibody Concentration:** 10.0 mg/ml (by UV absorbance at 280 nm,  $E^{0.1\%}$  of 1.4)

**Buffer:** 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

**Stabilizer:** None

**Preservative:** 0.01% (w/v) Sodium Azide

**Application(s):** Suitable for most immunological methods requiring purified immunoglobulins.

**Recommended Dilution(s):** Not specified. Optimal titers for applications should be determined by the researcher.

**Storage Conditions:** Store vial at 4° C prior to restoration. Restore with 5.0 ml of distilled water (or equivalent). Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.

**Purity:** This product was prepared from normal Rat serum by a multi-step process which includes delipidation and salt fractionation followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in precipitin arcs against anti-Rat Serum corresponding to gamma globulins.

**Immunogen:** Not applicable

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.

All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation.

All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 326, Gilbertsville, Pennsylvania, USA.