



ROCKLAND
Immunochemicals for Research

800-656-7625

fax 610-367-7825

Product Specification Sheet

Product: Cat Serum (lyophilized)

Code: D301-05

Lot # 13440

Size: 5.0 ml

Physical State: Lyophilized

Protein Concentration: 85 mg/ml (by Refractometry)

Buffer: None

Stabilizer: None

Preservative: 0.01% (w/v) Sodium Azide

Application(s): Suitable for use as a control or diluent in assays and most other immunological methods.

Recommended Dilution(s): Not specified. Optimal titers for applications should be determined by the researcher.

Storage Conditions: Store vial at 4° C prior to restoration. Restore with 5.0 ml of deionized water (or equivalent). Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.

Purity: This product was prepared from normal serum by delipidation and defibrination. Assay by immunoelectrophoresis resulted in multiple precipitin arcs against anti-Cat Serum.

Immunogen: Not applicable

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.