

## Product Specification Sheet

**Product:** Mouse Brain powder

**Code:** C115

**Lot #** 4005

**Size:** 1.0 g

**Physical State:** Powder

**Buffer:** None

**Stabilizer:** None

**Preservative:** None

**Preparation Method:** Tissue powders are produced by washing the various organs repeatedly in saline to remove blood followed by multiple precipitations in acetone until all lipid is removed. Acetone is removed by desiccation.

**Recommended Use:** Suitable for use as a specific adsorbent to remove unwanted reactivities or as a source material. Product provides a high nuclear content and is micronized to a fine powder.

**Storage Conditions:** Store vial at 2 to 8° C prior to opening. No expiration date is given for this product if properly stored.

**Strain:** Swiss-Webster

**Sex:** Mixed

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.