

## Product Specification Sheet

**Product:** Mouse Complement (lyophilized)

**Code:** C201-0005

**Lot #** 17992

**Size:** 5.0 ml

**Physical State:** Lyophilized

**Protein Concentration:** 85.0 mg/ml (by Refractometry)

**IgG Concentration:** not determined

**Buffer:** None

**Stabilizer:** None

**Preservative:** None

**Recommended Dilution(s):** Not specified.

**Storage Conditions:** Store vial at 4° C prior to restoration. Restore with 5.0 ml of diluent (supplied with product). Centrifuge product (at 2° to 8° C) if not completely clear. Aliquot contents and freeze at -70° C or below. **NOTE: COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY.** Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.

**Immunoelectrophoresis:** normal

**pH:** normal

**Hemoglobin:** normal

**Strain:** mixed breed

**Sex:** male

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.