

Product Specification Sheet

Product: Rabbit "A" (nesting) Complement (lyophilized)

Code: C202

Lot # 81992

Size: 5.0 ml

Physical State: Lyophilized

Protein Concentration: 85.0 mg/ml (by Refractometry)

IgG Concentration: not determined

Buffer: None

Stabilizer: None

Preservative: None

Recommended Dilution(s): Not specified.

Storage Conditions: Store vial at 4° C prior to restoration. Restore with 5.0 ml of diluent (supplied with product). Centrifuge product (at 2° to 8° C) if not completely clear. Aliquot contents and freeze at -70° C or below. **NOTE: COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY.** Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.

Immunoelectrophoresis: normal

pH: normal

Hemoglobin: normal

Strain: mixed breed

Sex: male

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.