

## Product Specification Sheet

**Product:** Guinea Pig Complement (fresh frozen)

**Code:** C300-0050

**Lot #** 32498

**Size(s):** 50 ml

**Physical State:** Fresh Frozen

**Protein Concentration:** 85.0 mg/ml (by Refractometry)

**IgG Concentration:** Not determined

**Buffer:** None

**Stabilizer:** None

**Preservative:** None

**Recommended Dilution(s):** Not applicable.

**Storage Conditions:** Store vial at -70° C prior to opening. **NOTE: COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY.** Avoid cycles of freezing and thawing.

**Immunoelectrophoresis:** Normal

**pH:** Normal

**Hemoglobin:** Normal

**Strain:** English Hartley short-haired

**Sex:** Male

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.