

Product Specification Sheet

Product: Rabbit "A" Nesting Complement (fresh frozen)

Code: C302

Lot # 65895

Size(s): 50 ml or 100 ml or 500 ml

Physical State: Fresh Frozen

Protein Concentration: 85.0 mg/ml (by Refractometry)

IgG Concentration: Not determined

Buffer: None

Stabilizer: None

Preservative: None

Recommended Dilution(s): Not specified.

Storage Conditions: Store vial at -70° C prior to opening. **NOTE: COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY.** Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of opening.

Immunoelectrophoresis: Normal

pH: Normal

Hemoglobin: Normal

Strain: Mixed breed

Sex: Mixed

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.