

Product Specification Sheet

Product: Bovine Plasma (non-sterile) in Citrate, Sodium

Code: D400-03-0050

Lot # 9033

Size: 50 ml

Physical State: Liquid

Protein Concentration: Normal

IgG Concentration: Not determined

Buffer: None

Stabilizer: None

Preservative: None

Recommended Dilution(s): Not specified.

Storage Conditions: Store vial at -20° C prior to opening. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of opening.

Sterilization: Not performed

Immunoelectrophoresis: Normal

pH: Normal

Hemoglobin: Normal

Strain: Mixed breed and sex

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.