

Product Specification Sheet

Product: Rabbit Gamma Globulin

Code: D610-0050

Lot # 13651

Size: 50 mg

Physical State: Lyophilized

Protein Concentration: 20 mg/ml (by UV absorbance at 280 nm, $E^{0.1\%}$ of 1.4)

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Stabilizer: None

Preservative: 0.01% (w/v) Sodium Azide

Application(s): Suitable for most immunological methods requiring purified immunoglobulins.

Recommended Dilution(s): Not specified. Optimal titers for applications should be determined by the researcher.

Storage Conditions: Store vial at 4° C prior to restoration. Restore with 2.5 ml of distilled water (or equivalent). Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.

Purity: This product was prepared from normal Rabbit serum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in precipitin arcs against anti-Rabbit Serum corresponding to gamma globulins.

Immunogen: Not applicable

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.