

## Product Specification Sheet

**Product:** DAB (3,3'-Diaminobenzidine) Substrate Tablets

**Code:** DAB-50

**Lot :** 12637

**Size :** 50 tablets (5 mg each)

**Physical State:** tablet

**Applications:** This product is designed to be used with the **MaxTag™** Immunoblotting kit. Dissolve one tablet in 10 ml of deionized water. Use liquid substrate within one hour. Immediately before use add hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) to 0.1% (v/v).

**Storage Conditions:** Store container at 2-8° C prior to opening. Keep dry. Warm to room temperature before opening. Expiration date is one (1) year from date of receipt.

**Shipping Conditions:** Use care. Ship as a hazardous material. DAB is a suspected carcinogen. Avoid exposure to eyes, skin and mucous membranes. Do not inhale or ingest. Wash hands thoroughly after using.

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.