

Product Specification Sheet

Product: Rabbit IgG control

Code: KIC003

Lot # 3872

Size: 0.5 ml

Physical State: Liquid (sterile filtered)

Antibody Concentration: 0.1 mg/ml (by UV absorbance at 280 nm)

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Stabilizer: 50% (v/v) Glycerol

Preservative: 0.01% (w/v) Sodium Azide, 0.01% (w/v) Gentamicin Sulfate

Application(s): suitable for use as a control or standard in most immunological assays.

Recommended Dilution(s): See Kit insert for complete instructions.

Storage Conditions: Store vial at 4° C. Dilute only prior to immediate use. Expiration date is one (1) year from date of opening.

Purity: This product was prepared from normal mouse serum by a multi-stage process which includes delipidation, salt fractionation and ion exchange chromatography. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit IgG and anti-Rabbit Serum.

Immunogen: not applicable

Note: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.