

## Product Specification Sheet

**Product:** Sterile Chicken Blood, Defibrinated

**Code:** R102-0100

**Lot #** 65987

**Size:** 100 ml

**Physical State:** Suspension

**Anticoagulant:** None

**Stabilizer:** None

**Preservative:** None

**Recommended Dilution(s):** Not specified.

**Storage Conditions:** Store product at 4° C. Use standard sterile techniques and sterile equipment at all times to maintain sterility.

**Sterilization:** This product was aseptically drawn using sterile equipment and standard sterile techniques. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.

**Expiration Date:** If properly stored and handled this product **may** be stable for up to one (1) month. Be advised that blood is a perishable product and exact shelf may depend on application.

**Animal(s):** normal healthy non-medicated adult Chicken of mixed sex and breed

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.