

## Product Specification Sheet

**Product:** Washed Pooled Guinea Pig Cells, 10%

**Code:** R402

**Lot #** 071598-G

**Size:** 50 ml or 100 ml

**Physical State:** Suspension

**Anticoagulant:** saline

**Stabilizer:** None

**Preservative:** None

**Recommended Dilution(s):** Not specified.

**Storage Conditions:** Store product at 4° C.

**Expiration Date:** If properly stored and handled this product **may** be stable for up to one (1) week. Be advised that blood is a perishable product and exact shelf may depend on application.

**Animal(s):** normal healthy non-medicated adult Guinea Pigs of mixed sex and breed

**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.